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| Method | GC/MS | | Total Results |
|----------------------|----------|----------|---------------|
| | Results | Positive | |
| MDMA Rapid Test | Positive | 20 | 22 |
| | Negative | 2 | 68 |
| Total Results | | 22 | 68 |
| % Agreement | | 90.9% | 97.1% |

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Non Cross-Reacting Compounds

| | | | |
|-------------------------|-------------------------|-------------------|------------------------|
| 4-Acetamidophenol | Dextromethorphan | Meprobamate | Procaine |
| Acetophenetidin | Diclofenac | Methamphetamine | Promazine |
| N-Acetylprocainamide | Diazepam | Methadone | Promethazine |
| Acetylsalicylic acid | Diffunisal | Methoxyphenamine | D,L-Propranolol |
| Aminopyrine | Digoxin | Methyphenidate | D-Propoxyphene |
| Amitypyline | Morphine | Morphine | D-Pseudoephedrine |
| Amobarbital | Diphenhydramine | 3-β-D-glucuronide | Quinacrine |
| Amoxicillin | 5,5 - Diphenylhydantoin | Morphine sulfate | Quinidine |
| Ampicillin | Doxylamine | Nalidixic acid | Quinine |
| I-Ascorbic acid | Egonine hydrochloride | Naloxone | Ranitidine |
| D-Amphetamine | Egoninemethylester | Naltrexone | Salicylic acid |
| D,L-Amphetamine sulfate | (-) -ψ-Ephedrine | Naproxen | Secobarbital |
| I-Amphetamine | [1R,2S](-) Ephedrine | Niacinamide | Serotonin |
| Apomorphine | l - Epinephrine | Nifedipine | (5-Hydroxytyramine) |
| Aspartame | Erythromycin | Nimesulidate | Sulfamethazine |
| Atropine | β-Estradiol | Norcocodol | Sulindac |
| Benzilic acid | Estrone-3-sulfate | Norethindrone | Sustiva |
| Benzolic acid | Ethyl-p-aminobenzoate | D-Norpropoxyphene | Temazepam |
| Benzoylcegonine | Fenoprofen | Noscapine | Tetracycline |
| Benzphetamine | Furosemide | D,l-Octopamine | Tetrahydrocortisone, |
| Bilirubin | Genticic acid | Oxalic acid | 3- Acetate |
| (±) - Brompheniramine | Hemoglobin | Oxazepam | Tetrahydrocortisone |
| Bupropion | Hydralazine | Oxolinic acid | 3-(β-D glucuronide) |
| Caffeine | Hydrochlorothiazide | Oxycodone | Tetrahydrozoline |
| Cannabidiol | Hydrocodone | Oxymetazoline | Thebaine |
| Cannabinol | Hydrocortisone | Papaverine | Theophylline |
| Chloralhydrate | O-Hydroxyhippuric acid | Penicillin-G | Thiamine |
| Chloramphenicol | p-Hydroxyamphetamine | Pentazocine | Trans-2- |
| Chloridiazepoxide | p-Hydroxy- | hydrochloride | phenylcyclopropylamine |
| Chlorothalidate | methamphetamine | Pentobarbital | Thioridazine |
| (±) - Chlorpheniramine | 3-Hydroxytyramine | Perphenazine | Tolbutamide |
| Chlorpromazine | Imipramine | Phencyclidine | Trazodone |
| Chlorquine | Iproniazid | Phenelzine | DJ-Tyrosine |
| Cholesterol | (±) - Isoproterenol | Phenobarbital | Triamterene |
| Clomipramine | Isosuprine | Phentemine | Trifluoperazine |
| Clonidine | Ketamine | Trans-2-phenyl | Trimethoprim |
| Cocaeethylene | Ketoprofen | cyclopropylamine | Trimipramine |

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| MDMA Concentration (ng/ml) | Percent Of Cut-off | n | Visual Result | |
|----------------------------|--------------------|----|---------------|----------|
| | | | Negative | Positive |
| 0 | 0 | 30 | 30 | 0 |
| 250 | -50% | 30 | 30 | 0 |
| 500 | Cut-off | 30 | 15 | 15 |
| 750 | +50% | 30 | 0 | 30 |
| 3X | 1500 | 30 | 0 | 30 |

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
| | |
|---------------------------------------|-----------|
| (±)3,4-Methylenedioxy-mpheta-mine HCl | 50 ng/ml |
| (±) 3,4-Methylenedioxyamphetamine HCl | 300 ng/ml |

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1. Tietz NW. Textbook of Clinical Chemistry. W.B. Saunders Company. 1986; 1735
2. Baselt RC. Disposition of Toxic Drugs and Chemicals in Man. 2nd Ed. Biomedical Publ., Davis, CA. 1982; 488

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