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Method	GC/MS		Total Results	
	Results	Positive		Negative
TCA Rapid Test	Positive	91	13	104
	Negative	5	141	146
Total Results		96	154	250
% Agreement		94.8%	91.6%	92.8%

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Nortriptyline Concentration (ng/mL)	Percent of Cut-off	n	Visual Result	
			Negative	Positive
0	0%	30	30	0
500	-50%	30	30	0
750	-25%	30	25	5
1,000	Cut-off	30	15	15
1,250	+25%	30	3	27
1,500	+50%	30	0	30
3,000	3X	30	0	30

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Compound	Concentration (ng/mL)	Compound	Concentration (ng/mL)
Nortriptyline	1,000	Imipramine	400
Nordoxepine	500	Clomipramine	50,000
Trimipramine	3,000	Doxepine	2,000
Amitriptyline	1,500	Maprotiline	2,000
Promazine	3,000	Promethazine	50,000
Desipramine	200	Perphenazine	50,000
Cyclobenzaprine	2,000	Dithiaden	10,000
Desipramine	200	Perphenazine	50,000
Cyclobenzaprine	2,000	Dithiaden	10,000

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WWW.PHARMEDCO.IR

Non Cross-Reacting Compounds			
Acetophenetidin	Dextromethorphan	Methadone	Phenylpropanolamine
N-Acetylprocainamide	Diazepam	D-methamphetamine	Prednisolone
Acetylsalicylic acid	Diclofenac	(L)-methamphetamine	Prednisone
Aminopyrine	Diffunisal	Methoxyphenamine	Procaine
Amobarbital	Digoxin	3,4-Methylenedioxyethyl- amphetamine	D,L-Propranolol
Amoxicillin	Diphenhydramine	amphetamines	D-Propoxyphene
Ampicillin	Doxylamine	(±) 3,4-Methylenedioxy- methamphetamine	D-Pseudoephedrine
L-Ascorbic acid	Ecgonine hydrochloride	methamphetamine	Quinidine
Apomorphine	Ecgonine methyl ester	Methylphenidate	Quinine
Aspartame	(1R,2S)-(-)-Ephedrine	Morphine-3-b-D- glucuronide	Ranitidine
Atropine	L-Ephedrine	Naloxone	Salicylic acid
D,L -Amphetamine	Erythromycin	Nalidixic acid	Secobarbital
L-Amphetamine	Ethyl-p-aminobenzoate	Naloxone	Serotonin
Benzilic acid	Fenfluramine	Naltrexone	(5-Hydroxytyramine)
Benzoic acid	Fenoprofen	Naproxen	Sulfamethazine
Benzoylcegonine	Furosemide	Niacinamide	Sulindac
Benzphetamine	Gentisic acid	Nifedipine	Temazepam
Bilirubin	Hemoglobin	Norcodeine	Tetracycline
(±)-Brompheniramine	Hydralazine	(-)-ψ- Ephedrine	Tetrahydrocortisone, 3 Acetate
Caffeine	Hydrochlorothiazide	Norethindrone	Tetrahydrocortisone 3 (b-D glucuronide)
Cannabidiol	Hydrocodone	D-Norpropoxyphene	Tetrahydrozoline
Cannabinol	Hydrocortisone	Noscapine	Thebaine
Chloralhydrate	p-Hydroxyamphetamine	D,L-Octopamine	Thiamine
Chloramphenicol	O-Hydroxyhippuric acid	Oxalic acid	Thioridazine
Chlordiazepoxide	3-Hydroxytyramine	b-Estradiol	Tolbutamide
Chlorothiazide	p-Hydroxy- methamphetamine	Oxycodone	Triamterene
(±) Chlorpheniramine	ibuprofen	Papaverine	Trifluoperazine
Chlorpromazine	(±)-Isoproterenol	Penicillin-G	Trimethoprim
Chlorquine	Isoxsuprine	Pentazocine	D, L-Tryptophan
Cholesterol	Ketamine	Pentobarbital	Tyramine
Clonidine	Ketoprofen	Phencyclidine	D, L-Tyrosine
Cocaine hydrochloride	Labetalol	Phenelzine	Uric acid
Codeine	Levorphanol	Phenobarbital	Verapamil
Cortisone	Loperamide	Phentermine	Oxazepam
(-) Cotinine	Meperidine	L-Phenylephrine	Zomepirac
Creatinine	Meprobamate	b-Phenylethylamine	
Deoxycorticosterone			

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1. Rose, J.B., Tricyclic antidepressants toxicity. J. Toxicity Clin. Toxicol. ۱۱, ۳۸۱-۴۰۲, ۱۹۷۷
2. Hawks RL, CN Chiang. Urine Testing for Drugs of Abuse. National Institute for Drug Abuse (NIDA), Research Monograph ۷۳, ۱۹۸۶